

Products for USP <797>

Compliance

From the leader in contamination control in critical environments



Williams Medical Company

A Division of Williams Medical Supply, Inc. Phone: 800.252.8646 - Fax: 714.961.8914 www.williamsmedical.com

Product Selection Guide

USP <797> MANDATE	RECOMMENDED TEXWIPE PRODUCTS		
	Sterile Environment	Non-sterile Environment	
Clean bottles, syringes, vials and other supplies transferred into the clean area	Pre-wetted Wipers Sterile AlphaSat [®] 10 Sterile TechniSat [®] Sterile PolySat [®]	Pre-wetted Wipers TechniSat [®] PolySat [®]	
	Dry Wipers and Sterile Alcohol SterileWipe [™] HS II SterileWipe [™] LP Sterile 70% Isopropanol	Dry Wipers and Sterile Alcohol QuanTex [™] AlphaWipe® Technicloth [®] VersaWipe® Sterile 70% Isopropanol	
Clean equipment: - Buffer room and anteroom carts - Walls and ceilings - Surfaces	Pre-wetted Wipers Sterile AlphaSat® 10 Sterile TechniSat® Dry Wipers and Sterile Alcohol	Pre-wetted Non-sterile Wipers TechniSat® PolySat® Dry Wipers and Sterile Alcohol	
- Floors - Shelving	SterileWipe [™] HS II SterileWipe [™] LP Sterile 70% Isopropanol	QuanTex [™] AlphaWipe® Technicloth [®] VersaWipe® Sterile 70% Isopropanol	
	Mops and Sterile Alcohol AlphaMop [™] AlphaMop [™] Covers Mini AlphaMop [™] ICT Sterile 70% Isopropanol	Mops and Sterile Alcohol AlphaMop [™] AlphaMop [™] Covers Mini AlphaMop [™] ICT Sterile 70% Isopropanol	
Documents in controlled environment	Cleanroom Notebooks TexWrite [®] Spirals TexWrite [®] Bound	Cleanroom Notebooks TexWrite [®] Spirals TexWrite [®] Bound	
Sanitize gloves	Pre-wetted Wipers Sterile AlphaSat® 10 Sterile TechniSat®	Pre-wetted Wipers TechniSat [®]	
	Dry Wipers and Sterile Alcohol SterileWipe [™] HS II SterileWipe [™] LP Sterile 70% Isopropanol	Dry Wipers Technicloth [®] VersaWipe®	
Clean floors, walls and ceilings Apply and remove disinfectants	Mops and Sterile Alcohol AlphaMop [™] Sterile 70% Isopropanol	Mops AlphaMop [™] AlphaMop [™] Covers AlphaMop [™] Pre-wetted Mop Covers	
Clean and disinfect: - Laminar flow hoods - Compounding Aseptic Isolators (CALs) - Biological Safety Cabinets (BSCs)	Mini AlphaMop [™] ICT CleanTips [®] Swabs Alpha® Belt Mopping System	Mini AlphaMop [™] ICT CleanTips® Swabs Alpha® Belt Mopping System	
Contamination control in buffer areas	Cleanroom Adhesive Mats FirstStep [®] Mat CleanStep [™] Mat CleanStep [™] Refills	Cleanroom Adhesive Mats FirstStep [®] Mat CleanStep [™] Mat CleanStep [™] Refills	
Environmental surface monitoring Cleaning validation Disinfection of ampoules and stoppers when wetted with IPA	CleanTips [®] Swabs	CleanTips [®] Swabs	



Sterile Dry Wipers

SterileWipe[™] HS II

Absorbent, hydroentangled cellulose/polyester wiper. ISO Class 6-7.

TX3210	12" x 12" wiper	100/bag	TX609	9" x 9" wiper	3
		20/pack 5 packs/bag	TX612	12" x 12" wiper	1

SterileWipe[™] LP

Knitted polyester wiper for use in critical sterile environments. ISO Class 3-5.

TX3211	9" x 9" wiper	100/bag 10/pack 10 packs/bag
TX3212	12" x 12" wiper	100/bag 20/pack 5 packs/bag



Dry Wipers

TechniCloth® Hydroentangled cellulose/polyester wiper, good strength and durability. ISO Class 6-7.

TX609	9" x 9" wiper	300/bag
TX612	12" x 12" wiper	150/bag

VersaWipe®

Hydroentangled cellulose/polyester wiper, high absorbency, best value. ISO Class 6-7.

TX622	12" x 12" wiper	150/bag
TX624	4" x 4" wiper	1,200/bag
TX629	9" x 9" wiper	300/bag

AlphaWipe®

Absorbent polyester wiper, low particles and fibers. ISO Class 5-6.

TX1004	4" x 4" wiper	300/bag
TX1009	9" x 9" wiper	150/bag
TX1013	12" x 12" wiper	75/bag

Vectra[®] QuanTex[™] Wipers

State of the art wiper provides maximum performance in critical environments. ISO Class 3-6.

TX1080 9" x 9" wiper 150/bag

Additional wiper products are available. Visit www.williamsmedical.com for more information.




Sterile Alcohol

Sterile 70% Isopropanol

Filtered to 0.2 microns and gamma irradiated.

TX3270	16 oz trigger-spray bottle	12/case
TX8270	8 oz. trigger-spray bottle	12/case
TX3290	1 gallon polybottle	4/case



Contamination Control Mats

FirstStep[®] Mat

Removes foot-borne particulates prior to entering the clean work area. An antimicrobial additive protects against microorganism growth. Adheres directly to floor.

AMA183613BB01	Blue	18" x 36" mat 60-layer	4 mats/box
AMA254513BB01	Blue	25" x 45" mat 60-layer	4 mats/box
AMA254513WW01	White	25" x 45" mat 60-layer	4 mats/box

Cleanroom Notebooks

TexWrite® 30 Spiral Notebooks

Economical choice for note taking and data recording. College-ruled pages.

TX5840	100 pages, 8.5" x 11"	10/box
TX5841	100 pages, 5.5" x 8.5"	10/box

TexWrite® 30 Bound Notebooks

Sewn binding prevents removal of pages without damage. Provides security of permanent, legally binding data records. College-ruled pages.

TX5807	Blue, 100 pages, 5.5" x 8.5"	10/box
TX5808	Blue, 100 pages, 8.5" x 11"	10/box

CleanStep[™] Mat

Provides the same protection properties as FirstStep against particulates and microorganisms with the ability and convenience to relocate the mat when necessary. Non-skid reusable frame.

AMQ263215WW31 W	 24" x 30" white refill mat	One frame and mat/box
	30-layer	

CleanStep[™] Mat Refills

Convenient replacement mats for use with the CleanStep frame.

AMS243017WW31	White	24" x 30" 30-layer	4 refills/box
AMT243013WW31	White	24" x 30" 60-layer	4 refills/box

Additional swabs are available. Visit www.williamsmedical.com for more information.



CleanTips[®] Swabs

Large Alpha[®] Swab

Knitted polyester head with polypropylene handle. Recommended for environmental surface monitoring and cleaning validation. Autoclavable. ISO Class 4-8.

TX714A

100/bag 10 bags/case

Total Swab Length: 5.01"

Small Alpha® Swab

Knitted polyester head with polypropylene handle. Wet with IPA for disinfecting stoppers on bags and vials before piercing with sterile needles and for necks of ampoules before they are broken. Autoclavable. ISO Class 4-8.

ТХ743В

500/bag 5 bags/case

Total Swab Length: 2.77"

Absorbond[®] Swab with Long Handle

Nonwoven polyester head with polypropylene handle. Recommended for surface sampling and cleaning hard-to-reach areas. May be wetted with IPA. Autoclavable. ISO Class 4-8.

TX762

100/bag 10 bags/case

Total Swab Length: 6.39"

Wrapped Cotton Swab

High sorption, 100% cotton fiber swab for routine cleaning and for spill clean-up. Autoclavable. ISO Class 6-8.

TX705

1,000/bag 5 bags/case

Total Swab Length: 5.98"

Additional swabs are available. Visit www.williamsmedical.com for more information.

Cleanroom Mops

AlphaMop™

Large polyester flat mop for wiping ceilings, walls and floors. ISO Class 4-8.

TX7108	Cleanroom mop with fiberglass handle and 15" x 7" single head assembly	1 mop/case
TX7118	Polyester replacement covers and foam pads for TX7108	150 covers/case 6 foam pads/case
TX718H	Honeycomb [®] 10 textured polyester covers for scrubbing	150/case

Mini AlphaMop[™] Isolator Cleaning Tool

Cleaning system for isolators, laminar flow hoods and biological safety cabinets.

TX7101	Isolator cleaning mop with two handles (18" and 24") (ISO Class 4-8)	1 mop/case
TX7111	TechniCloth [®] replacement covers (ISO Class 6-8)	150 covers/case 6 foam pads/case
TX7114	Polyester replacement covers (ISO Class 4-8)	150 covers/case 6 foam pads/case

Alpha[®] Belt Mopping System

Save time and money while eliminating buckets and carts in the controlled environment. ISO Class 4-8.

TX7211	Cleanroom belt, 16-oz. empty squirt bottle, two plastic rings, 20 disposable cleanroom bags	10/case	
TY7250	Alpha® Bolt Monning System	100/0000	

TX7250 Alpha[®] Belt Mopping System 100/case Refill Bags



Sterile Pre-Wetted Wipers with 70% IPA

Sterile TechniSat®

Hydroentangled cellulose/polyester wiper, pre-wetted with 70% IPA. ISO Class 6-7.

TX3214	9" x 11" wiper	50/flexpack
TX3217	9" x 11" wiper	20/flexpack

Sterile PolySat®

Polypropylene wipers pre-wetted with 70% IPA. ISO Class 6-7.

TX3213	9" x 11" wiper	50/flexpack
TX3216	9" x 11" wiper	20/flexpack

Sterile AlphaSat® 10

Polyester wiper with thermally sealed-borders, folded for small footprint, pre-wetted with 70% IPA. ISO Class 4-5.

TX3285

9" x 9" wiper

20/flexpack

Sterile AlphaSat[®] 10 with WFI Polyester wiper with thermally sealed-borders pre-wetted with 70% IPA/30% WFI. ISO Class 4-5.

TX3280

12" x 12" wiper

50/flatpack

Pre-Wetted Wipers with 70% IPA

PolySat[®]

Polypropylene wiper pre-wetted with 70% IPA, economical. ISO Class 6-7.

TX1051	9" x 11" wiper	50/flexpack
	•	•

TechniSat®

Hydroentangled cellulose/polyester wiper, pre-wetted with 70% IPA. ISO Class 6-7.

TX1041	9" x 11" wiper	70/flip-top tub
TX1045	6" x 8" wiper on roll	100/canister
TX1065	9" x 11" wiper	50/flexpack



ISO Class 3-5

Aseptic and Sterile Environments

Aseptic and sterile environments isolate the work area from the background environment.

ITW Texwipe has the tools, products and services to maintain the cleanliness of the work area.

TX7401 Mini Environment Cleaning and Maintenance Starter Kit

TX7101 — Mini AlphaMop[™] Isolator Cleaning Tool (ICT)

- Economical
- Easier to use and more ergonomic than wipers
 - · Cleans hard-to-reach areas quickly and easily
 - · Low, flat profile
 - Light weight can be used with one hand
 - Swivel joint flexible use
- Excellent for isolators, biosafety cabinets, glove boxes and laminar flow hoods
- Fits into most pass-throughs and ante-chambers
 - Easy to disassemble and reassemble
 - Two handles (18" and 24")
- Tolerates sterilization by autoclave, chemicals or gamma irradiation

TX3285 — Sterile AlphaSat® 10 with 70% IPA

- Sterile, individually folded, 9" x 9" wipers with a small package footprint (13" x 7"), saving valuable space
- Sealed-border wipers prevent fiber release
- Use to wipe down all surfaces, equipment and gloves
- Excellent for spill control in isolators, biosafety cabinets, glove boxes and laminar flow hoods
- Bag-within-a-Bag[®] cleanroom packaging
- Polyester wipers pre-wetted with 70% IPA
 - Safer and easier to use than a dry wiper and a bottle of IPA
 - · Eliminates handling of bottles of flammable liquids
 - Optimal concentration of IPA for killing microbes
- Resealable flexpack
- Use to remove cleaning and disinfectant residues



TX3211 — SterileWipe[™] LP Dry Wiper

- Sterilized for sterile and aseptic processes
- Ideal for cleaning, spill control and wiping gloves
- Bag-within-a-Bag[®] cleanroom packaging
- Good abrasion resistance
- Tolerates most solvents

This kit also includes the following literature:

- TechNote ICT-1: Cleaning of Compounding Aseptic Isolators (CAIs) and Class II Biological Safety Cabinets (BSCs)
- Brochures:
- Products for USP <797> Compliance CleanTips® Swabs
- MSDS: TX3285

For additional technical information go to www.texwipe.com



ISO Class 5

Controlled Environments

Controlled environments isolate the work area from the background environment.

ITW Texwipe has the tools, products and services to maintain the cleanliness of the work area.

TX7402 Mini Environment Cleaning and Maintenance Starter Kit

TX7101 —

Mini AlphaMop[™] Isolator Cleaning Tool (ICT)

- Economical
- Easier to use and more ergonomic than wipers
- Cleans hard-to-reach areas quickly and easily
- Low, flat profile
- Light weight can be used with one hand
- Swivel joint flexible use
- Excellent for isolators, biosafety cabinets, glove boxes and laminar flow hoods
- Fits into most pass-throughs and ante-chambers
 - Easy to disassemble and reassemble
 - Two handles (18" and 24")
- Tolerates sterilization by autoclave, chemicals or gamma irradiation

TX1084 — QuanSat[™] QuanTex[™] Wipers with 70% IPA

- Designed for ISO Class 3-7 critical environments
 - Wipe down all surfaces, equipment and gloves
 - Remove cleaning and disinfecting residues
 - · Excellent for spill control
- Resealable flatpack
- Bag-within-a-Bag[®] cleanroom packaging
- Sealed-border wipers prevent fiber release
- Polyester wipers pre-wetted with 70% IPA
 - · Safer and easier to use than a dry wiper and a bottle of IPA
 - Eliminates handling of bottles of flammable liquids
 - Optimal concentration of IPA for killing microbes

TX1009 — AlphaWipe® Dry Wiper

1084

Premium knit all-purpose heavy-weight polyester wiper

1009

AlphaWipe

- Designed for ISO Class 5 critical environments
- Cleanroom laundered and packaged
- Superior for spill control and general wiping
- Bag-within-a-Bag[®] cleanroom packaging
- Resists abrasion under rigorous use

This kit also includes the following literature:

- TechNote ICT-1: Cleaning of Compounding Aseptic Isolators (CAIs) and Class II Biological Safety Cabinets (BSCs)
- Product Brochures:
 - Products for USP <797> Compliance CleanTips® Swabs
- MSDS: TX1084

For additional technical information go to www.texwipe.com



A Division of Williams Medical Supply. In

ISO Class 7

Controlled Environments

Controlled environments isolate the work area from the background environment.

ITW Texwipe has the tools, products and services to maintain the cleanliness of the work area.

TX7403 Mini Environment Cleaning and Maintenance Starter Kit

TX7101 —

Mini AlphaMop[™] Isolator Cleaning Tool (ICT)

- Economical
- Easier to use and more ergonomic than wipers
- Cleans hard-to-reach areas quickly and easily
- Low, flat profile
- Light weight can be used with one hand
- Swivel joint flexible use
- Excellent for isolators, biosafety cabinets, glove boxes and laminar flow hoods
- Fits into most pass-throughs and ante-chambers
 - Easy to disassemble and reassemble
 - Two handles (18" and 24")
- Tolerates sterilization by autoclave, chemicals or gamma irradiation

TX1065 — TechniSat[®] with 70% IPA

- Low-linting wipers designed for cleanrooms
- Excellent for spill control and general wiping
- Cleaning the exterior of production and laboratory equipment
- Packaged in resealable flexpacks with small footprints (13" x 7")
- Polyester/cellulose wipers pre-wetted with 70% isopropyl alcohol (IPA)
 - Safer and easier to use than a dry wiper and a bottle of IPA
 - · Eliminates handling of bottles of flammable liquids
 - Optimal concentration of IPA for killing microbes

TX609 — TechniCloth® 9" x 9" Dry Wiper

- Cleanroom wiper
 - Better than cellulose wiper

⊺X"609

TechniCloth

- Low-linting
- · Highly absorbent
- Excellent wet and dry strength
- Excellent for spill control
- Fabricated from a blend of cellulose and polyester
 - · Highly absorbent properties of a natural fiber
 - Cleanliness and strength of a synthetic fiber
- Suitable for cleaning precision components and isolator apparatus

This kit also includes the following literature:

- TechNote ICT-1: Cleaning of Compounding Aseptic Isolators (CAIs) and Class II Biological Safety Cabinets (BSCs)
- Product Brochures:
- Products for USP <797> Compliance CleanTips® Swabs
- MSDS: TX1065



8

For additional technical information go to www.texwipe.com



Key products compounding pharmacies need in one convenient package!

USP <797> requires that the preparation of CSPs occur in ISO Class 5 (Class 100) environments, buffer areas of ISO Class 7 (Class 10,000), and ante-areas of ISO Class 8 (Class 100,000).

Achieving and maintaining cleanroom environments requires new procedures and protocols. ITW Texwipe[®] and Williams Medical have the key products compounding pharmacies need to meet these requirements.

ISO Class 7-8

Buffer Areas/ Ante-areas



the mop and laminated training posters are yours




WILLIAMS MEDICAL COMPANY Phone: 800-252-8646 Fax: 714-961-8914 www.williamsmedical.com **Unable Consistency.** Support.

www.texwipe.com

USP <797> Compliance Cleaning Kit

Kit contains 1 bag each of dry and pre-wetted wipers, a 16 ounce bottle of 70% isopropanol, 1 package of point-of-use pre-wetted wipers, 1 bag of mop covers and an adhesive floor mat with frame. The mop and training posters are included FREE!

тхкіт Complete kit





WILLIAMS MEDICAL COMPANY Phone: 800-252-8646 Fax: 714-961-8914 www.williamsmedical.com



Quality. Consistency. Support.

www.texwipe.com




Apparel and gloves for controlled environments



Williams Medical Company A Division of Williams Medical Supply, In

Tel 800.252.8646 • Fax 714.961.8914 • E-mail info@williamsmedical.com

Apparel and gloves for ISO 4/5 (Class 10/100) controlled environments



Have you also considered the following items?

Miscellaneous

- Pouches
- Sterilization Wrap
- Tubing
- Cuff Tape
- Sealers
- Cart/Equipment Covers

CardinalHealth

Gowning Step 3

Polypropylene

Hood

Gowning Step 5

Gowning Step 4 -Face Mask

Secure-Gard[®] Earloop

Gowning Step 6 -Coveralls Micro-Clean[®] Sterile

Gowning Step 1 -Scrubs (not visible) • Limited-Use Tops/Pants

Gowning Step 8 -

 Esteem[®] CP Sterile CP100BT[®] Sterile

Multi-Flex[®] Brand CC100[™]

Gloves

Step 2/7 -**Shoe Covers** Micro-Clean[®] Sterile Spunbonded Polypropylene

• SMS

and Masks


Bouffant (not visible)

Premium Bouffant Caps

• Comfort[™] Spunbonded

Micro-Clean[®] Sterile Combo

Micro-Clean[®] Sterile Hoods





Catalog Numbers ISO 4/5 (Class 10/100)

Gowning Step 1: Scrubs

Limited-Use Tops/Pants			
Cat. No.	Size	Qty. Cs	
Short Sleeve Top	os with V-Neck		
23502T	Medium	60/ea	
23503T	Large	60/ea	
23504T	X-large	60/ea	
23505T	XX-large	60/ea	
23506T	XXX-large	60/ea	
Long Sleeve Top	s with Round Neck		
23402LST	Medium	60/ea	
23403LST	Large	60/ea	
23404LST	X-large	60/ea	
23405LST	XX-large	60/ea	
23406LST	XXX-large	60/ea	
Pants with Elastic Waist			
23602PE	Medium	60/ea	
23603PE	Large	60/ea	
23604PE	X-large	60/ea	
23605PE	XX-large	60/ea	
23606PE	XXX-large	60/ea	

Gowning Step 2: Shoe Covers

Spunbonded Polypropylene			
Cat. No.	Size/Description	Qty. Cs	
2850	Universal	400/ea	
2852	Universal, skid-resistant	400/ea	
2854	X-large, skid-resistant	400/ea	
2859	XX-large, skid-resistant	400/ea	

Micro-Clean® Sterile 🔇

Cat. No.	Size	Qty. Cs
9823	Small	150/pr
9824	Medium	150/pr
9825	Large	150/pr
9826	X-large	150/pr
9827	XX-large	150/pr

SMS		
Cat. No.	Size/Description	Qty. Cs
Dura-Pak [™]		
4870	Universal	400/ea
4872A	Universal, skid-resistant	400/ea
4874	X-large, skid-resistant	400/ea
4875	XX-large, skid-resistant	400/ea



2850 Spunbonded Polypropylene Shoe Covers



9823 Micro-Clean® Sterile High-Top Shoe Covers



4870 Fluid-Resistant SMS Shoe Covers



Catalog Numbers ISO 4/5 (Class 10/100)

Gowning Step 3: Bouffant

Premium Bouffant Caps				
Cat. No.	Size	Color	Qty. Cs	
B21HD	21″ regular	Blue	1000/ea	
B24HD	24" large	Blue	750/ea	
W24HD	24" large	White	750/ea	
Comfort™ Spunbonded Polypropylene				
Cat. No.	Size	Color	Qty. Cs	
Dispenser Packed				
3244	21″ regular	Blue	1000/ea	

3244	21″ regular	Blue	1000/ea
3274	24″ large	Blue	750/ea
3277	25" X-large	Blue	750/ea
9305	21″ regular	White	1000/ea
9306	24" large	White	1000/ea
Bulk Packed			
3275	23″ medium	Blue	300/ea
9302	21″ regular	White	250/ea
9309	21″ regular	White	1000/ea
9303	24″ large	White	250/ea
9310	24″ large	White	1000/ea

Gowning Step 4: Face Mask

Micro-Clean® Sterile Hoods and Masks 🔇			
Cat. No.	Description	Qty. Cs	
9820	Hood with ties and mask	100/ea	
9810	Hood with ties and snap and mask	100/ea	
Secure-Gard® Earloop			
Secure-Gard [®] Ea Cat. No.	rloop Description	Qty. Cs	
	1	Qty. Cs 300/ea	
Cat. No.	Description		
Cat. No. AT7511	Description Blue earloop mask	300/ea	


Gowning Step 5: Hood

Micro-Clean [®]	Sterile Hoods and Masks 🔇	
Cat. No.	Description	Qty. Cs
9820	Hood with ties and mask	100/ea
9810	Hood with ties, snap and mask	100/ea

-FRI-

Gowning Step 6: Coveralls

Micro-Clean® Sterile 🖏				
Cat. No.	Size	Fits Chest	Fits Height	Qty. Cs
9831	Small	32″ to 36″	5´ to 5´7″	30/ea
9832	Medium	36″ to 38″	5′5″ to 5′10″	30/ea
9833	Large	38″ to 40″	5′8″ to 6′	30/ea
9834	X-large	40″ to 43″	5′11″ to 6′3″	25/ea
9835	XX-large	43" to 46"	6′1″ to 6′5″	25/ea
9836	XXX-large	46″ to 52″	6´2″ to 6´6″	29/ea
9837	4XL	52″ to 58″	6′3″ to 7′	29/ea
9838	5XL	58″ to 64″	6′5″ to 7′5″	29/ea





3244 Comfort[™] Spunbonded Polypropylene Bouffant



9820 Micro-Clean® Sterile Hood and Mask



AT7511 Secure-Gard® Earloop Mask



9831 Micro-Clean® Sterile Coveralls



A Division of Williams Medical Supply, In

Catalog Numbers ISO 4/5 (Class 10/100)

Gowning Step 7: Sterile Shoe Covers

Micro-Clean® Sterile 🖏			
Cat. No.	Size	Qty. Cs	
9823	Small	150/pr	
9824	Medium	150/pr	
9825	Large	150/pr	
9826	X-large	150/pr	
9827	XX-large	150/pr	

Gowning Step 8: Gloves (Double)

Esteem® CP 😻 🙍			
Cat. No.	Size	Qty. Cs	
2Y72PL55	51/2	200/pr	
2Y72PL60	6	200/pr	
2Y72PL65	6 ¹ / ₂	200/pr	
2Y72PL70	7	200/pr	
2Y72PL75	7 ¹ / ₂	200/pr	
2Y72PL80	8	200/pr	
2Y72PL85	8 ¹ / ₂	200/pr	
2Y72PL90	9	200/pr	

CP100BT® Latex 🔇

Cat. No.	Size	Qty. Cs
2Y5000	51/2	200/pr
2Y5001	6	200/pr
2Y5002	6 ¹ / ₂	200/pr
2Y5003	7	200/pr
2Y5004	7 ¹ / ₂	200/pr
2Y5005	8	200/pr
2Y5006	8 ¹ / ₂	200/pr
2Y5007	9	200/pr

Multi-Flex[®] Brand CC100[™] Latex

Cat. No.	Size	Qty. Cs
2Y1600	5 ¹ / ₂	200/pr
2Y1601	6	200/pr
2Y1602	6 ¹ / ₂	200/pr
2Y1603	7	200/pr
2Y1604	7 ¹ / ₂	200/pr
2Y1605	8	200/pr
2Y1606	8 ¹ / ₂	200/pr
2Y1607	9	200/pr



2Y72PL55 Esteem® CP Sterile Polyisoprene Clean Process Glove



2Y5000 CP100BT® Sterile Latex Clean Process Glove



2Y1600 Multi-Flex[®] Brand CC100[™] Latex Cleanroom Glove



Apparel and gloves for ISO 6 (Class 1000) controlled environments

For Sterile or Nonsterile



Have you also considered the following items?

Hood

- Comfort[™] Hood
- Tyvek[®] Hood
- Micro-Clean® Sterile

Coveralls

SMS

High-Top Shoe Covers

- Micro-Clean® 2-1-2
- Micro-Clean[®] Sterile

Scrubs

• Limited-Use Tops/Pants

Miscellaneous

- Pouches
- Sterilization Wrap
- Tubing
- Cuff Tape
- Sealers
- Cart/Equipment Covers



Gowning Step 2

Polypropylene

Gowning Step 3 – Face Mask

Secure-Gard[®] Earloop

Secure-Gard[®] N95

Gowning Step 4

Material with Snaps

Gowning Step 5 -

Micro-Clean[®] Sterile

Gowning Step 7

Gowning Step 1/6 Shoe Covers • Dura-Trac™ • SMS

Micro-Clean[®] Sterile

Sleeves

Gloves • Duraprene[™] CP • CP100[™] Nitrile Ambi

Tyvek[®] Sleeves

Tyvek[®] Sterile Premier-Style

Lab Coat/Frock

SMS

Premium Bouffant Caps

• Comfort[™] Spunbonded

Bouffant



Catalog Numbers ISO 6 (Class 1000)

Gowning Step 1: Shoe Covers

Dura-Trac™		
Cat. No.	Size	Qty. Cs
9326	Universal	300/ea
9327	X-large	300/ea
SMS		
Cat. No.	Size/Description	Qty. Cs
Dura-Pak™		
4870	Universal	400/ea
4872A	Universal, skid-resistant	400/ea
	· · · · , · · · · · ·	
4874	X-large, skid-resistant	400/ea

Micro-Clean® Sterile 🔇

Cat. No.	Size	Qty. Cs
9823	Small	150/pr
9824	Medium	150/pr
9825	Large	150/pr
9826	X-large	150/pr
9827	XX-large	150/pr

Gowning Step 2: Bouffant

Premium Bouffant Caps

Cat. No.	Size	Color	Qty. Cs
B21HD	21" regular	Blue	1000/ea
B24HD	24" large	Blue	750/ea
W24HD	24" large	White	750/ea

Comfort[™] Spunbonded Polypropylene

		-	
Cat. No.	Size	Color	Qty. Cs
Dispenser Packee	d		
3244	21" regular	Blue	1000/ea
3274	24" large	Blue	750/ea
3277	25" X-large	Blue	750/ea
9305	21" regular	White	1000/ea
9306	24" large	White	1000/ea
Bulk Packed			
3275	23″ medium	Blue	300/ea
9302	21" regular	White	250/ea
9309	21" regular	White	1000/ea
9303	24" large	White	250/ea
9310	24" large	White	1000/ea



9326 Dura-Trac[™] Shoe Covers



4870 Fluid-Resistant SMS Shoe Covers



9823 Micro-Clean® Sterile High-Top Shoe Covers



B21HD Premium Bouffant Cap



3244 Comfort[™] Spunbonded Polypropylene Bouffant



Catalog Numbers ISO 6 (Class 1000)

Gowning Step 3: Face Mask

Secure-Gard® Earloop					
Cat. No.	Cat. No. Description Qty. C				
AT7511	AT7511 Blue earloop mask				
AT7511-G	Green earloop mask	300/ea			
AT7511-PK	Pink earloop mask	300/ea			
AT7511-W	White earloop mask	300/ea			
N95 Respirator					
Cat. No.	Description	Qty. Cs			
N95-S	Small BX20	120/ea			
N95-ML MD/LG BX20		120,00			
N95-ML	MD/LG BX20	120/ea			
N95-ML FRN95-S	MD/LG BX20 2 Straps, SM BX20				
		120/ea			
FRN95-S	2 Straps, SM BX20	120/ea 120/ea			
FRN95-S FRN95-ML	2 Straps, SM BX20 2 Straps, MD/LG BX20	120/ea 120/ea 120/ea			

Gowning Step 4: Lab Coat/Frock

SMS			
Cat. No.	Size	Length	Qty. Cs
2401F	Medium	41″	100/ea
2402F	Large	42″	100/ea
2403F	X-large	44″	100/ea
2404F	XX-large	46″	100/ea
2405F	XXX-large	48″	100/ea
2406F	4X-large	50″	100/ea

Tyvek® Sterile Premier-Style Frocks with Snaps 🕉

Cat. No.	Size	Fits Chest	Fits Height	Qty. Cs
A6085-31S	Medium	32" to 36"	5' to 5'7"	25/ea
A6085-32S	Large	36" to 38"	5′5″ to 5′10″	25/ea
A6085-33S	X-large	38" to 40"	5'8" to 6'1"	25/ea
A6085-34S	XX-large	40" to 43"	5'11" to 6'3"	25/ea
A6085-35S	XXX-large	43" to 46"	6'1" to 6'5"	25/ea
A6085-36S	4X-large	46" to 52"	6'2" to 6'6"	25/ea



AT7511 Secure-Gard® Earloop Mask



N95 Respirator



2401F Fluid-Resistant Frock



A6085-31S Tyvek[®] Sterile Premier-Style Frock with Snaps



Catalog Numbers ISO 6 (Class 1000)

Gowning Step 5: Sleeves

Tyvek [®] Sleeves	Š	
Cat. No.	Size	Qty. Cs
9402	18″ L	60/ea
9403N	18″ L	300/ea
Micro-Clean [®] S	terile	
Cat. No.	Size	Qty. Cs
A6128-26	18″ L	50/pr

Gowning Step 6: Sterile Shoe Covers

Micro-Clean® Sterile 🔇		
Cat. No.	Size	Qty. Cs
9823	Small	150/pr
9824	Medium	150/pr
9825	Large	150/pr
9826	X-large	150/pr
9827	XX-large	150/pr

Gowning Step 7: Gloves

Duraprene [™] C	Duraprene™ CP 🖏 🧟			
Cat. No.	Size	Qty. Cs		
2Y72PN55	5 ¹ / ₂	200/pr		
2Y72PN60	6	200/pr		
2Y72PN65	6 ¹ / ₂	200/pr		
2Y72PN70	7	200/pr		
2Y72PN75	7 ¹ / ₂	200/pr		
2Y72PN80	8	200/pr		
2Y72PN85	8 ¹ / ₂	200/pr		
2Y72PN90	9	200/pr		

CP100[™] Nitrile Ambi 逆

Cat. No.	Size	Qty. Cs
2Y1849	X-small	1000/ea
2Y1850	Small	1000/e
2Y1851	Medium	1000/ea
2Y1852	Large	1000/ea
2Y1853	X-large	1000/ea



9402 Tyvek[®] Sleeves



A6128-26 Micro-Clean® Brand Tyvek® Sterile Sleeve



9823 Micro-Clean® Sterile High-Top Shoe Covers



89027-156 Duraprene[™] CP Latex-Free Glove



2Y1849 CP100[™] Nitrile Ambi Cleanroom Glove



Apparel and gloves for ISO 7/8 (Class 10,000/100,000) controlled environments

For Sterile or Nonsterile

Gowning Step 2 Bouffant

- Premium Bouffant Caps
- Comfort[™] Spunbonded Polypropylene

Gowning Step 3 Lab Coat/Frock

- SMS
- Tyvek[®] Sterile Premier-Style Material With Snaps

Gowning Step 4 – Nonsterile Gloves

- Esteem[®] Stretchy Nitrile
- Multi-Flex® Clean Process
 Ambi Latex
- CP100[™] Nitrile
- Duraprene[™] CP Sterile

Gowning Step 1

- Shoe Covers
- Dura-Trac™
- SMS



Have you also considered the following items?

Sleeves

- Tyvek[®] Sleeves
- Micro-Clean[®] Sterile

Hood

- Comfort[™] Hood
- Micro-Clean[®] Sterile

Coverall

• SMS

Face Mask

- Secure-Gard[®] Earloop
- Secure-Gard[®] N95

High-Top Shoe Covers

- Micro-Clean[®] Sterile
- Micro-Clean®

Scrubs

 Limited-Use Tops/Pants

Miscellaneous

- Pouches
- Sterilization Wrap
- Tubing
- Cuff Tape
- Sealers
- Cart/Equipment
- Covers





Catalog Numbers ISO 7/8 (Class 10,000/100,000)

Gowning Step 1: Shoe Covers

Dura-Trac™		
Cat. No.	Size	Qty. Cs
9326	Universal	300/ea
9327	X-large	300/ea
SMS		
Cat. No.	Size/Description	Qty. Cs
Dura-Fit [™] (Dispens	ser Box)	
4850	Universal	400/ea
4852	Universal, skid-resistant	400/ea
4854	X-large, skid-resistant	400/ea
4855	XX-large, skid-resistant	400/ea
Dura-Pak™		
4870	Universal	400/ea
4872A	Universal, skid-resistant	400/ea
4874	X-large, skid-resistant	400/ea
4875	XX-large, skid-resistant	400/ea

Gowning Step 2: Bouffant

Premium Bouffant Caps

Cat. No.	Size	Color	Qty. Cs
B21HD	21" regular	Blue	1000/ea
B24HD	24" large	Blue	750/ea
W24HD	24" large	White	750/ea

Comfort[™] Spunbonded Polypropylene

connore spansonaca rospropyiene						
Cat. No.	Size	Color	Qty. Cs			
Dispenser Packed	Dispenser Packed					
3244	21" regular	Blue	1000/ea			
3274	24" large	Blue	750/ea			
3277	25" X-large	Blue	750/ea			
9305	21" regular	White	1000/ea			
9306	24" large	White	1000/ea			
Bulk Packed						
3275	23″ medium	Blue	300/ea			
9302	21" regular	White	250/ea			
9309	21" regular	White	1000/ea			
9303	24" large	White	250/ea			
9310	24" large	White	1000/ea			



9326 Dura-Trac[™] Shoe Covers



4850 Fluid-Resistant SMS Shoe Covers



B21HD Premium Bouffant Cap



3244 Comfort[™] Spunbonded Polypropylene Bouffant



Catalog Numbers ISO 7/8 (Class 10,000/100,000)

Gowning Step 3: Lab Coat/Frock

SMS			
Cat. No.	Size	Length	Qty. Cs
2401F	Medium	41″	100/ea
2402F	Large	42″	100/ea
2403F	X-large	44″	100/ea
2404F	XX-large	46″	100/ea
2405F	XXX-large	48″	100/ea
2406F	4X-large	50″	100/ea

Tyvek® Sterile Premier-Style Frocks with Snaps 🕉

Cat. No.	Size	Fits Chest	Fits Height	Qty. Cs
A6085-31S	Medium	32" to 36"	5' to 5'7"	25/ea
A6085-32S	Large	36" to 38"	5′5″ to 5′10″	25/ea
A6085-33S	X-large	38" to 40"	5′8″ to 6′1″	25/ea
A6085-34S	XX-large	40" to 43"	5'11" to 6'3"	25/ea
A6085-35S	XXX-large	43" to 46"	6′1″ to 6′5″	25/ea
A6085-36S	4X-large	46" to 52"	6'2" to 6'6"	25/ea



2401F Fluid-Resistant Frock



A6085-31S Tyvek[®] Sterile Premier-Style Frock with Snaps



Catalog Numbers ISO 7/8 (Class 10,000/100,000)

Gowning Step 4: Gloves

Duraprene [™] CP	• 🔇 💮	
Cat. No.	Size	Qty. Cs
2Y72PN55	5 ¹ / ₂	200/pr
2Y72PN60	6	200/pr
2Y72PN65	6 ¹ / ₂	200/pr
2Y72PN70	7	200/pr
2Y72PN75	7 ¹ / ₂	200/pr
2Y72PN80	8	200/pr
2Y72PN85	8 ¹ / ₂	200/pr
2Y72PN90	9	200/pr

CP100[™] Nitrile Ambi 👾

Cat. No.	Size	Qty. Cs
2Y1849	X-small	1000/ea
2Y1850	Small	1000/ea
2Y1851	Medium	1000/ea
2Y1852	Large	1000/ea
2Y1853	X-large	1000/ea

Esteem[®] Stretchy Nitrile

Cat. No.	Size	Qty. Cs
2Y1870T	X-small	1000/ea
2Y1871T	Small	1000/ea
2Y1872T	Medium	1000/ea
2Y1873T	Large	1000/ea
2Y1874T	X-large	1000/ea

Multi-Flex® Clean Process Ambi Latex

Cat. No.	Size	Qty. Cs
2Y1700	Small	1000/ea
2Y1701	Medium	1000/ea
2Y1702	Large	1000/ea



89027-156 Duraprene™ CP Latex-Free Glove



2Y1849 CP100[™] Nitrile Ambi Cleanroom Glove



2Y1870T Nuthera Glove



2Y1700 Multi-Flex® Clean Process Ambi Latex Glove



Compounded Steri ow-Risk Leve Preparations

Contemporation and a state EMPTY VIAL

Compounded Sterile



Compounded Sterile High Risk Level ²reparations



Low Risk

products (CSP's). All employees compounding low risk preparations must have successfully passed the low risk media-fill test procedure at least preparation of low risk compounded sterile The Low Risk Media Test Kit is designed to document correct personal technique and The Low Risk Media Test Kit consists of effective environmental control for the annually.

3 - 50mL Sterile Empty Vials 3 sets of 4 - 5mL Test Media

instructions and the following:

24



RUX HIGH

ELT MERIA





TUNER RUNER





Medium Risk

passed the medium risk media-fill test procedure medium risk preparations must have successfully preparation of medium risk compounded sterile The Medium Risk Media Test Kit is designed to products (CSP's). All employees compounding The Medium Risk Media Test Kit consists of document correct personal technique and effective environmental control for the at least annually or more frequently. nstructions and the following:



High Risk

addition, pass a media-fill test that represents a products (CSP's). All employees compounding high risk preparations must pass the low risk preparation of high risk compounded sterile The High Risk Media Test Kit is designed to document correct personal technique and media-fill test procedures annually and in The High Risk Media Test Kit consists of high-risk level every six months by each effective environmental control for the instructions and the following: authorized person.

9 - 10mL Sterile Empty Vials 1 - 100cc Test Medium





www.williamsmedical.com Tel 800.252.8646 • Fax 714.961.8914 • E-mail info@williamsmedical.com

VCA INTRAVENOUS QUALITY ASSURANCE

A Division of Williams Medical Supply, In Williams Medical Company

About Williams Medical Company

Established in 1978 Williams Medical Company distributes quality pharmaceuticals and pharmacy disposable to a variety of healthcare professionals and medical institutions. We are a stocking distributor who maintains a physical inventory for every product line that we offer.



Williams Medical Company is an authorized distributor for every manufacturer that we represent.



All of the products that we offer are purchased directly from the manufacture. We meet and exceed each manufacturer's stringent requirements with respect to accountability, handling, storage, and shipment. Our temperature-controlled warehouses are at a constant temperature throughout the year. We retain "lot number histories" for every product sold. Sales reports are provided to each manufacture in order to insure that all of their respective sales representatives are informed about the sales and to insure that they receive their appropriate credit. Further, we are a licensed pharmacy wholesaler for every State in the Nation.

Experts in Sales and Service, we maintain a professional staff of full time salespeople and associates. People who are at your door routinely and consistently to attend to your every need. Order the wrong item? A product going out of date? An Emergency need? We are there to address and resolve those issues for you.

Need Same Day Service? Try Us! Any delivery can be made in Southern California within three hours. Any delivery can be made in the nation in 12-24 hours. That's 7 Days a week, 24 hours a day for your convenience!

Order Online from 40 different manufacturers-24 Hours a Day! Think of the time that you'll save by only placing one order for all of the products that you need from 40 different manufactures. Forty less Purchase Orders! Forty less orders to receive and process. Forty less invoices for your accounting department to pay, when one payment will do. While Online, you can review the latest FDA recalls, obtain CE credits, review your order history, and see that's new!

Our Commitment to you is to provide a level of service and reliability that you will grow to expect. We will save you time ordering and receiving goods. We will save you money accounting for these orders. We will save you the waste of having to destroy expired pharmaceuticals by addressing these issues before they go out of date. We will offer products at a fair price. We pledge that you will always speak with someone who is sincerely dedicated to meeting your needs and expectations.



MEETING YOUR USP 797 NEEDS



www.williamsmedical.com - Phone: 800.252.8646 - Fax: 714.961.8914 - A Division of Williams Medical Supply, Inc.